

The Test Culture of some Microorganisms with Whale Meat Peptone

Shichiro Akiya and Fumiharu Kobo

This experiment was carried out to determine the difference of value of utilization between the samples and market-peptones. Market peptones were used as follow;

Teruuchi's, Witte's, Carneau's and Chaptoteaut's Peptone.

Experimental

Incubation-condition ;

Culture media ; Peptone-water prepared with the some kinds of peptones as above described.

Test-strains ; B. coli communis, B. typhi, Pneumo coccus, B. dysenterie Minoda, V. cholerae, B. subtilie M. B. 25, Staph. coccus aureus Terajima, B. pyocyanus.

Incubation ; 37°C., 24 hrs.

Result ;

The amount of increased bacillus, the formation of dyestuff and H₂S are shown as follow :

B. Coli communis

Peptone	Hours	Indol-formation				
		Ehrlich's Method		Kitazato's Method		
	6	12	18	24	Order	
Sample	+	+	#	#	1	# 2
Teruuchi	#	#	#	#	2	# 1
Carneau	+	#	#	#	3	+ 5
Witte	±	+	+	#	5	# 4
Chaptoteaut	+	#	#	#	4	# 3

B. typhi

Peptone	Hours	Indol-formation				
		6	12	18	24	Order
Sample	+	#	#	#	2	
Teruuchi	+	#	#	#	1	
Carneau	+	#	#	#	4	
Witte	+	+	+	#	5	
Chaptoteaut	+	#	#	#	3	

Pneumococcus 1

Peptone	Hours	Indol-formation				
		6	12	18	24	Order
Sample	#	#	#	#	4	
Teruuchi	#	#	#	#	1	
Carneau	#	#	#	#	5	
Witte	#	#	#	#	3	
Chaptoteaut	+	#	#	#	2	

B. dysenteriae Minoda

Peptone	Hours	6	12	18	24	Order
Sample		+	+	+	#	4
Teruuchi		#	#	#	#	1
Carneau		+	#	#	#	2
Witte		±	+	+	#	5
Chaptoteaut		+	#	#	#	3

V. cholerae O.

Peptone	Hours	6	12	18	24	Order	Indol-formation	Order
Sample		#	#	#	#	1	+	5
Teruuchi		+	+	#	#	4	+	4
Carneau		#	#	#	#	3	+	5
Witte		+	+	#	#	5	+	2
Chaptoteaut		±	+	+	#	2	#	1

The experiments above described led summarily to the results as follow :

Peptone	Order of Bac.-growth					Order of Indol-formation		
Sample	1	1	2	4	4	2	3	5
Teruuchi	1	1	1	2	4	1	3	4
Carneau	2	3	3	3	3	2	5	5
Witte	2	5	5	5	5	2	4	4
Chaptoteaut	2	2	3	3	4	1	1	3

(i) The amount of increased Bac. Incubation ; 24 hrs, Peptone-water, pH 7.0

Strain	B. subtilis M. B. 15 B. pyocyanus				Staph. aureus Terajima	
Peptone		Order		Order		Order
Teruuchi	#	1	#	1	#	1
Witte	#	4	#	2	#	1
Sample	#	2	#	1	#	1
Carneau	#	3	#	1	#	1
Chaptoteat	#	3	#	1	#	1

(ii) *B. typhi*; H₂S-Formation tested by Lead-acetate-paper. pH 6.8

(iii) Formation of dyestuff and Indol. B. Pyocyanus. Incubation; 36hrs, pH 7.0

Hours	24	36	Reaction	Dyestuff	Indol-formation	Growth
Peptone			Peptone		Ehrlich's Method	Kitazato's Method
Teruuchi	#	#	Teruuchi	-	#	#
Witte	+	#	Witte	#	+	#
Sample	-	-	Sample	-	#	#
Carneau	#	#	Carneau	#	#	#
Chaptoteaut	-	-	Chaptoteaut	+	+	#

Authors finally deeply thank Prof. Shogo Hosoya, Prof. Michizo Asano and Mr. Goro Urakubo of Infections Disease Institute of Tokyo University.

(Pharmaceutical department, the Faculty of Medicine, University of Tokyo and Whales Research Institute.)